Automation of Quantitative NMR Spectroscopy of Polyolefins in Industry

G. Hubner, M. Parkinson

InnoTech Operational Support, Borealis Polyolefine GmbH, Linz, Austria

H. Beiler, S. Steuernagel

Bruker Biospin GmbH, Rheinstetten, Germany

Pacifichem 2015 - Session: NMR Spectroscopy of Polymers Honolulu, USA 19.12.2015



Keep Discovering

Outline

- Polyolefin chain microstructure and quantitative NMR spectroscopy
- Automation of polyolefin solution-state NMR spectroscopy
 spectral analysis
- Automation of polyolefin melt-state NMR spectroscopy
 - sample exchange
- Summary and Conclusion

Polyolefin Chain Microstructure & Quantitative NMR Spectroscopy

Polyolefin Chain Microstructure & NMR



Polyolefin Chain Microstructure

- Isotacticity
- Comonomer content
- Short-chain branch distribution (kCbb, kCtotal)
- Regio-defect content (2,1e, 2,1t, 3,1 ...)
- Tacticity distribution (dyad, triad, pentad)
- Unsaturation content
- Unsaturation distribution
- Comonomer sequence distribution (dyad, triad)
- End-group distribution
- Long-chain branch content
- End-group content (Mn)
- Grafting content
- Oxidation











All possible from a single ¹³C NMR spectrum with no calibration (primary method)



Solution-State NMR Spectroscopy of Polyolefins





Polyolefins: 10 mm probehead and tube, TCE-d₂, 120°C, ¹³C detection



Modern Automated NMR spectroscopy

- Modern NMR spectrometers are very complicated but easy to use
 - Seriously amazing amount of automation (very high degree of abstraction)
 - Key aim of vendors: Better ease of use leads to more installations
- 1. Sample Track / LIMS interface
- 2. IconNMR / Automation interface
- 3. TopSpin / GUI interface
- 4. Backend / API
- 5. Embedded spectrometer control computer
- 6. Advanced timing hardware
- 7. Signal generation > Amplifies > Pulses > Probehead > Receiver > FID

automated by macros (linear CLI) automated by AU API (full C-program)





- Direct relative microstructure quantification from singles 1D spectrum
 - Quantitative NMR spectroscopy is easy in theory
 - Quantitative ¹³C{¹H} NMR not easy in practice
 - Quantitative ¹³C{¹H} NMR of polyolefins is difficult (valid)
- Must consider whole process of quantitative NMR





Quantitative NMR Spectroscopy Process



- Method development process works in reverse
 - fix what do we really want to know; define problem
 - how do we calculate it from a truly quantitative spectrum; implement quant method
 - asses if spectrum is quantitative by varying conditions and adjust



Automated Solution-State NMR Spectroscopy of Polyolefins

Automated Polyolefin Solution-State NMR



Custom Spectral Analysis: quant AU Programs

- Bruker TopSpin command line AU: quant-pb
- 1. Pass command line flags for flow control (21e noxx noends)
- 2. Integrate
- 3. Assign representative integrals/do compensation
- 4. Calculates primary microstructure quantities (contents, distributions)
- 5. Calculated secondary microstructure quantities (statistics)
- 6. Stochastic fitting (linear, non-linear)
- 7. Report in various formats (screen, PDF, CSV, TXT, LIMS)
- Multiple methods (literature regions/equations, custom in-house)
 - quant-ep-83 currently determines ethylene content in 26 different ways per spectrum
- Approximately 2500 line of code in C
 - readability more important than optimisation



Reporting: quant-eh-32 PDF



Topspin > processpe; quant-eb xx xex; pdf quant; catcsv



Reporting: quant-eh-32

##[Comonomer Content]####################################											
	.00	mol%	0.25	RNDL 0.29 0.3	ASTM 0.26 0.3	1 5 wt%	0.73 0.7	RNDL 0.86 0.9	ASTM 0.78 0.8		comonomer content
##[Sho	rt Chain Bra	anching]######	#####	*#########	##########	*######	#####[()2]##	
			Cbb	CHn	CH2			RNDL	ASTM		SCB content
	.00	butyl	1.23	1.22	1.23	3 mol%C6 ->	brB4/Cbb	1.45	1.30		
	.0	branches	1.2	1.2	1.2			1.4	1.3		
	.00	methyl	0.00!	0.00!	0.00)!					
	.0	branches	0.0!	0.0!	0.0!						
	0.0	total	1 23	1 22	1 23	3					
	.0	branches	1.2	1.2	1.2	sum butyl	and methy	vl branc	hes		
						_					
	.00	methyl aroups	4.25 4.3	4.23 4.2	4.26 4.3	all per de	efined 100	0 carbo	n atoms		
		920420									



epo	rting:	quant	t -eh-3	2				
Comonom	er Distribu	tion]####	##########	#########	##########	######	####[03]##	
				•••••				
Rand	lall89a tria	d method fo	or systems	with > 3	mol% copc	lymer		triad comonomer
								distribution
two	ranges C1 &	C2 (better	r)	single range C (worse)				distribution
EEE	98.95 %	XXX 0	.00 %	EEE	98.96 %	XXX	0.00 %	
EEX	0.76 %	EXX 0	.00 %	EEX	0.76 %	EXX	-0.01 %	
XEX	0.00 %	EXE 0	.29 %	XEX	0.00 %	EXE	0.29 %	
E	99.71 %	*X* 0.	.29 %	*E*	99.71 %	*X*	0.29 %	
prop	erties deri	ved from ti	riad distr	ibution				
		Σζ.	F		V	F		triad derived
	longth	A 1 0 267				Е ГО О		
sey.	number	1.0 20.			1.0 20	0 1		properties
seq.	number	0.5	J.4 1 0		0.3	0.4		
run		0.0	J . O	_	21 0	0.0		
CIUS	iter iux	0.0			21.9			
hran	ch defects	and realat	tive comon	omer dist	ribution w	ia inte	rnal method	
	len, dereets							
	branches		defec	ts	comono	mer		defect distribution
	br kCbb	0	kCbb	00	9			
EXE	Bu 1.2	100.0	1.2	100.0	100.0	1		
EXX	Bu 0.0!	0.0!	0.0!	0.0!	0.0	!		
XEX	Bu 0.0!	0.0!	0.0!	0.0!	0.0	!		
XXX	Bu 0.0!	0.0!	0.0!	0.0!	0.0	!		
В1	Me 0.0!	0.0!	0.0!	0.0!				
ALL	1.2!	100.0	1.2!	100.0	100.0)		

BOREALIS

Reporting: quant-eh-32

			Cbb	CHn	CH2	split					end groups
	.00	sat	3.02 3.0	3.01 3.0	3.02 3.0	95.8	% sat	:vin =	23.0		
	.00	unsat	0.13	0.13	0.13	4.2	% vin	:sat =	0.0		
	.0		0.1	0.1	0.1						
##	[Average	Chain	Propertie	s]#####	#####	#########	########	########	#####[05]##	
		cha	in length	634.2	carbo	ns ap	prox. Mn	8.92	k (8922.2)		average chain
	degree of bran	polym nches	erisation per chain	317.1	monom branc	ers hes	Mn MRL4	8.92 354.5	k (8922.2)		properties
			_								
##	E Signal-	to-Noi	se Ratios]#######	#####	#########	######### 	*#######	#####[06]##	
	si	te r	ange		SNR	%SD					SNR &
	(*1	dd+ 3 84 3	0.9 - 29.1 8 7 - 37 9	687 1	3.7 2 4	0.015		XX XEX	flag = 0		runtime flags
	al	B4 3	5.0 - 34.5	2	2.7	4.404		XXX	flag = 0		randing hage
	ć	aa 4	0.9 - 40.0		1.6	62.381		B1	flag = 0		
	11	вв 2 В1 2	4.9 - 24.4 0.4 - 19.9		1.2	84.325 61.391	use	useAlpha AllAlpha	fiag = 0 flag = 0		
- 6											







Automated Melt-State NMR Spectroscopy of Polyolefins

Automated Polyolefin Melt-State NMR



Melt-state MAS Vs Solution-State NMR Spectra



- Non standard technique based on solid-state NMR spectroscopy method
 - Compromise between sensitivity (higher) and resolution (lower)
 - Key technique is Magic-Angle Spinning (MAS)

M. Parkinson, et al. ACS Symp. 2011, 1077, 24, 401



Melt-state Magic-Angle Spinning (MAS) stator





- Standard: solid-state NMR
 - Place bulk sample in 7 mm zirconia ceramic rotor
 - Seal using a ceramic cap with find (turbine)
 - Levitate on high-pressure gas bearing inside stator
 - Drive rotation to 4 kHz with high-pressure gas
- Non-standard: melt-state NMR
 - Melt-the sample 150-180°C
 - Keep sample in rotor while under MAS
 - complex, high risk = expensive (9 kEUR)

z 57.74° Υ θm

G. Hatfield, et al. Anal. Chem. 1995, 67, 3082
K. Klimke, 2006 Ph.D. Thesis, University of Mainz, Germany
M. Pollard, et al. Macromolecules 2004, 37, 813
K. Klimke, et al. Macromol. Chem. Phys. 2006, 207, 382
M. Parkinson, et al. ACS Symp. 2011, 1077, 24, 401



Quantitative Melt-State NMR Spectroscopy Process



- In the world of solid-state NMR spectroscopy
 - 7 mm rotor packing is difficult (symmetric weight distribution) but in melt no issue!
 - MAS spin-up automated but not stable for melt-state (crash)
 - sample exchange very rarely automated (long experiments, different setup needed)





Major Challenges of Automating Melt-State MAS NMR

Difficulty

- Automated melt-state sample exchange: roadblock!
 - Very difficult: 2 years
- Automatic melt-state MAS Spin-up: roadblock!
 - Medium difficulty: 6 months
- Automatic tune/match (ATM) MAS probehead
 Medium difficulty: 3 days (validated compromise)
- Automatic spectral analysis (Don't measure what you can't analysis)
 - Easy as already done! 1 day

All solution also applicable to true solid-state NMR spectroscopy!





Automated Melt-State Sample Exchange

- Is there a current commercial solution available from Bruker Biospin?
 - 180 position sample changer for 7 mm MAS rotors capable of handling 180°C rotors
 - No: and was too complicated for in-house or automation specialist (ask vendor!)
- Can Bruker Biospin custom build hardware development
 - No: No market, non-standard setup, high-risk, little return on invested time/money (business)
 - However, did have large capacity changer for 4 mm MAS rotors designed for < 0°C!</p>
 - Significant modification needed
 - True collaboration with Bruker co-hardware development with open knowledge sharing
 - Bruker agreed to modify hardware for 7 mm (no liability, very special SAT)
 - Borealis agreed to handle all hot/heavy rotor modifications and melt-MAS spin-up
- Why did vendor proceed?
 - We had the technical competence
 - Realistic expectations: we thought it might not work either!



Bruker Custom 7 mm SamplePro: Trays



NOTE: rotors are not pre-heated and not store hot, just returned hot after measurement



Bruker Custom 7 mm SamplePro: Overview



- In theory is simple, we just blow the rotor down plastic tube
- In practice very difficult to get to work



Summary & Conclusions

Summary & Conclusion



Automated Polyolefin NMR Spectroscopy: Positive

- Automated spectral analysis for quantitative polyolefin NMR (solution & melt-state)
 - Higher quality
 - Very high repeatability (most recent melt-state method had a CV 0.1%!)
 - Removal of user dependent results (no more manual phasing/baseline/integral bias)
 - Full version control of programs QM/QA compliant (continuous improvement)
 - Rapid method development via automated analysis (1-2 weeks with LIMS)
 - Directly comparable method (same logic) e.g. EB, EH, EBH
- Automated melt-state NMR
 - Higher throughput (20+ to 160 quantitative ¹³C spectra/week @ same quality)
 - Faster delivery of project goals, new opportunities, all bench scale samples
 - all high-throughput Multi Modal Reactor (MMR) samples via ¹³C (8x20 min = 3 hr / day)
 - Higher efficiency/hardware utilisation (60% to 96%)
 - NMR spectrometers need to run 24/7/365 to break-even



Automated Polyolefin NMR Spectroscopy: Negative

- Time: all automation was done in parallel to our daily work!
- Higher overhead!
 - More time is spent keeping devices working now combined with less time to think
- Rapid backlog generation
 - When things break-down sample stack up very quickly
- Forget the complexity
 - Requestors: get accustomed to high-quality rapid feedback unrealistic expectations regarding method development
 - Laboratory: over time you forget how your automation works (comment the code!) miss small things (need checklists)
 - Specialists: underestimate work between proof-of-concept and implementation!





Summary

- Implemented advanced automation for both solution-state and non-standard melt-state NMR spectroscopy of polyolefins
- Implemented modular system for automated polyolefin NMR spectral analysis
- Successfully automated non-standard melt-state MAS NMR technique
 reached potential through full automation (long journey)
- Only possible with strong open vendor interaction
 - Custom hardware sample exchanger co-development (Bruker BioSpin IP)
 - Custom software melt-state MAS spin-up process (Borealis IP)



Acknowledgments

- InnoTech Operational Support (IOS) Management
 - Jochen Berrens, Jens Reussner

IOS NMR group

... the money & faith

- Bruker Biospin Solid-State NMR/automation
 - Steffan Steuernagel, Heiko Beiler, Jens Abromeit Michael Engelhardt, Rudolf Schöps

Gerhard Hubner, Lisa Steiner, Elena Pomakhina



... the trust

... the hard work

... the samples



- Polyolefin R&D groups
 - SPPR-PP, SPPR-PE, SPPR-HPPE, POLY, CAT-D, CAT-P

Thank you

A project by Borealis AG. The ideas documented in this presentation are the sole property of Borealis AG, and are subject to current copyright laws. Unauthorised use, reproduction in whole or in part, as well as transmission to third parties is not permitted.



NMR Automation Humor

- Q: What's worse than asking your management for money to automate a device?
- A: Asking for four times that amount one year later just for sample containers to run the device at capacity!



7 mm zirconia MAS rotor/cap set = 1 kEUR each



Future Perspectives: Standardisation

- Standardisation of polyolefin quantification by ¹³C NMR
 - currently no robust or universal standards (technically complex)
 - disconcerting when you consider how important tacticity of iPP is!
- Standardisation is in the interest of both industry and academia:
 - for Industry it would allows direct comparison for IPR/patent issues
 - for academia ensures repeatability/trust within academia and industry
- Technically and politically complex
 - need to standardise whole process on expensive heterogeneous equipment i.e. sample prep > measurement > processing> analysis > reporting
 - Sharing of proprietary knowledge/know-how from industry is complicated

